



CuraPro

Partial foot relief, post-operative and therapeutic shoes

**Orthopaedic aid
directory Germany:**

Forefoot relief shoes
31.03.03.5011

Heel relief shoes
31.03.03.5012

Longterm post-operative shoes
31.03.03.4001/2

CLINICALLY
TESTED



THANNER
SCHÄFTE
und mehr

Quality of Life

Movement boosts your quality of life and spells out independence. Unfortunately, a wide variety of injuries, illnesses and fractures can limit your mobility.

The solution

CuraPro - the longterm post-operative shoes, partial foot relief shoes and therapeutic shoes from THANNER enable you to stay in motion again.

Approved medical product

THANNER CuraPro special shoes are a class I MDD article 2(a) approved medical product, with CE labelling. THANNER also had their shoes tested by TÜV Rhineland [German Association for Technical Inspection] so that they are deemed an approved medical device. The material has been tested for safety, carcinogenicity and toxicology.

Made in Germany

All CuraPro special shoes are made by our 205 employees in our own facility in Höchstädt on the Danube, Germany. That's German quality.

Clinically tested

The effectiveness of the CuraPro longterm post-operative shoes and partial foot relief shoes have been proven by several clinical studies.

Orthopaedic aid number Germany

partial foot relief shoes:

Cabrio/Campus/Cambridge/Canterbury: 31.03.03.5011
Colmar/Comet/Concept: 31.03.03.5012

longterm post-operative shoes:

Calypso/Colombo: 31.03.03.4001
Caprice/Clavigo: 31.03.03.4002

Other post-operative shoes and partial foot relief shoes:
Certification procedures underway

17 years of experience

For more than 17 years, THANNER has been manufacturing post-operative shoes and partial foot relief shoes here. The experience we have gathered and close relationship we have built up to prominent medical specialists and orthopaedic shoe technicians has spelled out continuous improvements. Today you can choose from genuine state-of-the-art products.

Dealers

CuraPro special shoes can only be purchased through our dealers (orthopaedic shoe technology, qualified medical supplies retailer, orthopaedic technology). The orthopaedic shoemaker or orthopaedic technician will guarantee optimum orthopaedic aid fit and the option of modifying the insert template or shoe.

For left and right feet

CuraPro special shoes are the only ones produced for left and right feet. That means that they have optimum hold so that wounds can heal more quickly.

Tailor-made


You can always order longer velcro closings, other upper girth, different shoesizes and other special productions for a small extra charge.

Quality and environmental management

THANNER is certified in conformity with DIN EN ISO 13485:2010-01 for medical products and DIN EN ISO 14001 for the environment. This guarantees you highest-quality approved materials and structured production processes. Sustainability in product manufacture and environmental needs are important to us.

Pressure measurement

We use the Pedar pressure measurement system from Novel for manufacturing or coming up with advancements for CuraPro special shoes. This enables us to create maximum pressure relief for standard shoes.



Introduction


01

Partial foot relief shoes
Ladies and Men

03-06

Longterm
post-operative shoes
Ladies and Men

07-08

Pentagon insert templates
Ladies and Men

09-10

Heilbronn rocker shoes
Ladies and Men

11-12

Wound net shoes
Ladies and Men

13

Therapeutic shoes
Ladies and Men

14

Special and
prosthetic socks
Ladies and Men

15-17

So light ...

CuraPro by THANNER



Partial foot relief shoes

concentrate the step load on the strongest part of the foot. Impact protection has the purpose of protecting the toes and the cushion insert in the forefoot can be milled through for total relief. Beyond this, there is wide opening and individual lock regulation; reliable shoe fastening on the foot's step surface. Hypoallergenic materials tested for pollutants or biocompatibility. One of the few shoes that is produced for left and right (not single-balled) which provides better heel hold.

Left or right as single shoes available.

Special characteristics

- Partial foot relief shoes are produced for left and right feet (not single-balled). This is the reason why these shoes give better hold and wounds or injuries can heal more quickly.
- class I MDD article 2(a) medical product with CE labelling
- cushion insert 1 cm thick for removing
- can be supplied piecewise
- pollution-tested materials, some also with biocompatibility test

Custom-made productions

We can manufacture custom-made productions at any time in a wide variety of designs (extra charge).

- other colours and materials
 - different width or size for left and right
 - over or undersizes
 - orthotic compensation in the shoe or on the outsoles
- Custom-made productions are excluded from exchange or return.

Indication

Post-operative or other wound processes, also treated with extensive (even cushioned) bandages that heal over a longer period of time and allow at least partial foot relief when walking – even outside.

Prescription

- 1 piece of forefoot relief shoes
CuraPro (DIABETICUS+RHEUMATICUS®)
Cabrio/Campus/Cambridge/Canterbury
- 1 piece of heel relief shoes
CuraPro (DIABETICUS+RHEUMATICUS®)
Colmar/Comet/Concept

Orthopaedic aid number Germany

Cabrio/Campus/Cambridge/Canterbury: 31.03.03.5011

Colmar/Comet/Concept: 31.03.03.5012

Please use the insole template listed here to find the correct size and width!

| | | | | | |
|-------------------------|----------------|----------------|---------------|--|---|
| Shoe last 900 | 34 - 42 | 41 - 47 | LADIES MEN | COLMAR, COMET, CONCEPT | Insole outline Heel relief shoe, Longterm post-operative shoe |
| | SS, S, M | | | | |
| Shoe last 910 | 34 - 42 | 41 - 47 | LADIES MEN | CABRIO, CAMPUS CAMBRIDGE, CANTERBURY | Insole outline Forefoot relief shoe, Prosthetic shoe |



Heel relief shoes

IN THE ORTHOPAEDIC AID AND MEDICINE CATALOGUE GERMANY:

Colmar/Comet/Concept: 31.03.03.5012



You get better fit and quicker wound healing
by producing shoes for **left and right**.

| | |
|----------|---|
| 1 | <p>COMET 900</p> <p>Heel relief shoe 4.2010 • Width M • left 4.2011 • Width M • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin thermoplastic shoe counter</p> |
|----------|---|

| | |
|----------|--|
| 2 | <p>CONCEPT 900</p> <p>Heel relief shoe 4.2037 • Width M • left 4.2038 • Width M • right Size 34/35 - 46/47 (Double sizes) marine/ozean Space/Aquarin thermoplastic shoe counter</p> |
|----------|--|

| | |
|----------|---|
| 3 | <p>COLMAR 900</p> <p>Heel relief shoe 4.2012 • Width M • left 4.2013 • Width M • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin Lining Diasoft-3D Stiffening around the ankle</p> |
|----------|---|

Hallux valgus foot relief shoe



1

Forefoot relief shoes

IN THE ORTHOPAEDIC AID AND MEDICINE CATALOGUE GERMANY:

Cabrio/Campus/Cambridge/Canterbury: 31.03.03.5011



2



3



4

RAIN CAP

5

6



7

You get better fit and quicker wound healing by producing shoes for **left and right**.

| | |
|---|--|
| 1 | CAMBRIDGE 910 Hallux valgus foot relief shoe 4.2003 • left 4.2004 • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin |
|---|--|

| | |
|---|---|
| 3 | CAMPUS 910 Forefoot relief shoe 4.2005 • left 4.2006 • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin Lining Diasoft-3D |
|---|---|

| | |
|---|--|
| | RAIN CAP for forefoot relief shoe Size 34/35 - 46/47 (Double sizes) |
| 5 | 4.9859 • left 4.9860 • right oliv Space |
| 6 | 4.9862 • left 4.9863 • right marine Space |

| | |
|---|---|
| 2 | CABRIO 910 Forefoot relief shoe 4.2000 • left 4.2001 • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin Lining pocket on the heel |
|---|---|

| | |
|---|--|
| 4 | CANTERBURY 910 Forefoot relief shoe 4.2039 • left 4.2040 • right Size 34/35 - 46/47 (Double sizes) marine/ozean Space/Aquarin |
|---|--|

| | |
|---|--|
| 7 | CANTERBURY 910 Forefoot relief shoe 4.2041 • left 4.2042 • right Size 34/35 - 46/47 (Double sizes) oliv/schwarz Space/Hydro-Nappa |
|---|--|



Longterm post-operative shoes

(used over several months)

are used to encourage early ambulatory capability with extensive bandages needed on feet and legs combined with the rest of the foot being resilient and for protecting bandages from outside influences. Beyond this, there is wide opening and individual lock regulation, upper material adjustment and bottom adjustment option. The materials that are used are safe in terms of carcinogenicity and toxicology, hypoallergenic material = washable.

Left or right as single shoes available.

Special characteristics

- Post-operative shoes produced for left and right feet (not single-balled). This is why these shoes have better hold and wounds and injuries can heal more quickly.
- class I MDD article 2(a) medical product with CE labelling
- cushion insert 1 cm thick for removing
- can be supplied piecewise
- pollution-tested materials, some also with biocompatibility test

Custom-made productions

We can manufacture custom-made productions at any time in a wide variety of designs (extra charge).

- other colours and materials
- different width or size for left and right
- over or undersizes
- shortening compensation in the shoe or on the outsoles
- sole reinforcements made of carbon (with amputations) or fibreglass

Custom-made productions are excluded from exchange or return.

Indication

Post-operative or other wound processes, also treated with extensive (even cushioned) bandages that heal over a longer period of time and allow at least partial foot relief when walking - even outside.

Prescription

1 piece of longterm post-operative shoes
CuraPro (DIABETICUS+RHEUMATICUS®)

Orthopaedic aid number Germany

Calypso/Colombo: 31.03.03.4001

Caprice/Clavigo: 31.03.03.4002

Please use the insole template listed here to find the correct size and width!

| | | | | | |
|-------------------------|----------------|----------------|---------------|--|--|
| Shoe last 900 | 34 - 42 | 41 - 47 | LADIES MEN | CALYPSO, CAPRICE, COLOMBO, CLAVIGO | Insole template Heel relief shoe, Longterm post-operative shoe |
| | SS, S, M | | | | |
| Shoe last 910 | 34 - 42 | 41 - 47 | LADIES MEN | CABRIO, CAMPUS CAMBRIDGE, CANTERBURY | Insole template Forefoot relief shoe, Prosthetic shoe |

Longterm post-operative shoes

You get better fit and quicker wound healing by producing shoes for **left and right**.



IN THE ORTHOPAEDIC AID AND MEDICINE CATALOGUE GERMANY:

Calypso/Colombo Nr. 31.03.03.4001
Caprice/Clavigo Nr. 31.03.03.4002

Prosthetic shoe

INSOLE LENGTH


| SIZE | CHOPART | LISFRANC |
|-------|---------|----------|
| 34/35 | 14.0 cm | 16.9 cm |
| 36/37 | 14.8 cm | 17.3 cm |
| 38/39 | 15.7 cm | 18.6 cm |
| 40/41 | 16.6 cm | 19.4 cm |
| 42/43 | 17.4 cm | 20.3 cm |
| 44/45 | 18.3 cm | 21.2 cm |
| 46/47 | 19.1 cm | 22.0 cm |



| | |
|----------|--|
| 1 | COLOMBO 900 Longterm post-operative shoe 4.2045 • left 4.2046 • right Size 34-47 • Width M marine/ozean Space/Aquarin thermoplastic shoe counter |
| 2 | CALYPSO 900 Longterm post-operative shoe 4.1020 • left 4.1022 • right Size 34-47 • Width M ozean Aquarin thermoplastic shoe counter |

| | |
|----------|--|
| 3 | CAPRICE 900 Longterm post-operative shoe 4.1030 • left 4.1032 • right Size 34-47 • Width M ozean Aquarin Lining pocket on the heel |
| 4 | REGENKAPPE for longterm post-operative shoes 4.9864 • left 4.9865 • right Size 34/35 - 46/47 (Double sizes) marine Space |

| | |
|----------|---|
| 5 | CLAVIGO 900 Longterm post-operative shoe 4.2043 • left 4.2044 • right Size 34-47 • Width M marine/ozean Space/Aquarin |
| 6 | CHOPART 910 Prosthetic shoe 4.2015 • left 4.2016 • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin Lining pocket on the heel |
| 7 | LISFRANC 910 Prosthetic shoe 4.2017 • left 4.2018 • right Size 34/35 - 46/47 (Double sizes) ozean Aquarin Lining pocket on the heel |



The Pentagon insert template can be ordered as accessories for the Canterbury (forefoot relief shoe), Clavigo and Colombo (longterm post-operative shoes) shoe models.

You can alternate the inserts to lower Shore hardness 18° to relieve pressure.

You also have the option of completely removing the element for even greater relief.

Designed for a maximum body weight of 100 kg.

Fast pressure relief without any tools.

A great tip from Thanni:

Pressure relief without any tools!




There are 5 alternating elements in Shore 30° and Shore 18° for equipping the Pentagon insert template series



"Pentagon" Insert templates

| | |
|----------|--|
| 1 | PENTAGON 910 FOR CANTERBURY 4.4016 • left 4.4017 • right Size 34/35 - 46/47 (Double sizes) 3 mm Lunairmed Shore 18° 10 mm Dialux Shore 30° 10 mm Dialux Shore 18° |
|----------|--|

| | |
|----------|--|
| 2 | PENTAGON 900 FOR CLAVIGO/COLOMBO 4.4018 • left 4.4019 • right Size 34 - 47 3 mm Lunairmed Shore 18° 10 mm Dialux Shore 30° 10 mm Dialux Shore 18° |
|----------|--|



Heilbronn rocker shoe (Longterm post-operative shoes on a biomechanical basis)

are used to encourage early ambulatory capability with extensive bandages needed on feet and legs combined with the rest of the foot being resilient and for protecting bandages from outside influences. Beyond this, there is wide opening and individual lock regulation, upper material adjustment and outsole adjustment option.

The materials that are used are safe in terms of carcinogenicity and toxicology, hypoallergenic material = washable (biocompatibility test).

Left or right as single shoes available.

Special characteristics

- Heilbronn rocker shoes are produced for left and right feet (not single-balled). This is why these shoes have better hold and wounds or injuries can heal more quickly.
- class I MDD article 2(a) medical product with CE labelling
- can be supplied piecewise
- pollution-tested materials, some also with biocompatibility test
- can be worn with bandaging material
- minimises risk of thrombo-embolism

Custom-made productions

We can manufacture custom-made productions at any time in a wide variety of designs (extra charge).

- other colours and materials
- different width or size for left and right
- over or undersizes
- shortening compensation in the shoe or on the outsoles
- sole reinforcements made of carbon (with amputations) or fibreglass
- "Weitungen" on the lasts (more measure than width WW, Hallux valgus or claw toes)

Custom-made productions are excluded from exchange or return.

Indications

Toe arthrosis plastic (especially with Hallux valgus and hammer toes), corrective osteotomies of the toes and metatarsal bones, AVK Stad. IV, diabetic foot (osteolyses and ulcerations), post-operative phase after amputation of toes and forefoot, fractures of the toes, stable-exercise toe and metatarsal osteosyntheses, operations on the soft part of the forefoot

Prescription

1 piece of longterm post-operative shoes (Heilbronn rocker shoe)
CuraPro (DIABETICUS+RHEUMATICUS®)

Please use the insole template listed here to find the correct size and width!

| | | | | | | |
|-------------------------|----------------|----------------|---------------|------------|---|---|
| Shoe last 990 | 35 - 42 | 41 - 48 | LADIES | MEN | BREISGAU, HEIDELBERG, HEILBRONN, MARKGRÖNINGEN | Insole template Heilbronn rocker shoe, Wound net shoe |
|-------------------------|----------------|----------------|---------------|------------|---|---|



You get better fit and quicker wound healing by producing shoes for **left and right**.



| | |
|----------|--|
| 1 | HEIDELBERG 990 4.2027 • left 4.2028 • right Size 35/36 - 47/48 (Double sizes) schwarz Haftflex (Biocompatibility test) thermoplastic shoe counter |
|----------|--|

| | |
|----------|---|
| 2 | HEILBRONN 990 4.2020 • left 4.2021 • right Size 35/36 - 47/48 (Double sizes) schwarz Haftflex (Biocompatibility test) thermoplastic shoe counter |
|----------|---|

| | |
|----------|--|
| 3 | BREISGAU 990 4.2025 • left 4.2026 • right Size 35/36 - 47/48 (Double sizes) schwarz Vectra (Biocompatibility test) thermoplastic shoe counter |
|----------|--|

| | |
|----------|--|
| 4 | MARKGRÖNINGEN 990 4.2035 • left 4.2036 • right Size 35/36 - 45/46 (Double sizes) oliv/schwarz Space/Hydro-Micrella thermoplastic shoe counter |
|----------|--|

| | |
|----------|--|
| 5 | POLYETHYLEN - INSERT 4.4007 Size 35/36 - 47/48 (Double sizes) |
|----------|--|

| | |
|----------|--|
| 6 | POLYETHYLEN - INSERT WITH NEOPRENE COVER 4.4006 Size 35/36 - 47/48 (Double sizes) |
|----------|--|

Developed in cooperation with Physician Dr. P. Bause, Department of Trauma and Reconstructive Surgery, SLK-Kliniken Heilbronn GmbH



Wound net shoe

(longterm post-operative shoes used over several months)

You get better fit and quicker wound healing by producing shoes for **left and right**.

Indications

Post-operative or other wound processes, also treated with extensive (even cushioned) bandages that heal over a longer period of time and allow at least partial foot relief when walking - even outside.

Left or right as single shoes available.

Prescription

1 piece of longterm post-operative shoes
CuraPro (DIABETICUS+RHEUMATICUS®)
(Wound net shoe)



| | |
|----------|---|
| 1 | WUNDNETZ "M" 900 (WOUND NET SHOE) 4.2030 • left 4.2031 • right Size 35/36 - 47/48 (Double sizes) schwarz Haftflex (Biocompatibility test) thermoplastic shoe counter |
|----------|---|

| | |
|----------|---|
| 2 | WUNDNETZ "W" 990 (WOUND NET SHOE) 4.2032 • left 4.2033 • right Size 35/36 - 47/48 (Double sizes) schwarz Haftflex (Biocompatibility test) thermoplastic shoe counter |
|----------|---|

Please use the insole template listed here to find the correct size and width!

| | | | | | |
|-------------------------|----------------|----------------|---------------|--------------------------------------|--|
| Shoe last 900 | 34 - 42 | 41 - 48 | LADIES MEN | WUNDNETZ "M" (WOUND NET SHOE) | Insole template Heel relief shoe, Longterm post-operative shoe |
| | SS, S, M | | | | |
| Shoe last 990 | 35 - 42 | 41 - 48 | LADIES MEN | WUNDNETZ "W" (WOUND NET SHOE) | Insole template Heilbronn rocker shoe, Wound net shoe |

Therapeutic shoes



You get better fit and quicker wound healing by producing shoes for **left and right**.



Velcro closings can be removed

As a preventative measure with long positioning and ulcer positioning. The Haftflex material can be cut and washed.

| | |
|----------|--|
| 1 | CAMPANILE 950 4.1637 • left 4.1638 • right Size 34-47 oliv Space Sole: Race schwarz |
|----------|--|

| | |
|----------|--|
| 2 | CAMPINO 950 4.1635 • left 4.1636 • right Size 34-47 schwarz Haftflex (Biocompatibility test) thermoplastic shoe counter Sole: Race schwarz |
|----------|--|

| | |
|----------|---|
| 3 | BETTLAGERUNGS-SCHUH (DEKUBITUS) 4.1615 Size 34-38, 39-42, 43-47 single-ball (piecewise available) schwarz Haftflex with 20 mm cushion |
|----------|---|

Please use the insole template listed here to find the correct size and width!

| | | | | | | |
|-------------------------|----------------|----------------|---------------|------------|------------------------------------|--|
| Shoe last 950 | 34 - 42 | 41 - 47 | LADIES | MEN | CAMPANILE CAMPINO | Insole template Therapeutic last |
|-------------------------|----------------|----------------|---------------|------------|------------------------------------|--|

Group of persons

Special socks for ladies and men

Colours

Available in 5 different colours:
white, beige, marine, grey and black

Shoe sizes

- 35 to 38
- 39 to 41
- 42 to 43
- 44 to 45
- 46 to 47

Material composition

- 80 % cotton
- 20 % polyamide

Special characteristics

Available in 2 designs

- a) thick with terry
- b) thin

- machine washable at 60 °C, rinse carefully without fabric softener
- class I MDD article 2(a) medical product with CE labelling
- clinically tested

Custom-made socks productions

We can manufacture custom-made productions at any time in a wide variety of designs (extra charge).

- with individual logo or without any logo
- produced turned inside out (seams to the outside)
custom-made productions are excluded from exchange or return
- customised according to measurements





Special socks

| SIZES | 1 weiß | 2 beige | 3 marine | 4 grau | 5 schwarz |
|----------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| 35 - 38 | 5.3000 thick 5.3100 thin | 5.3010 thick 5.3110 thin | 5.3030 thick 5.3130 thin | 5.3040 thick 5.3140 thin | 5.3020 thick 5.3120 thin |
| 39 - 41 | 5.3001 thick 5.3101 thin | 5.3011 thick 5.3111 thin | 5.3031 thick 5.3131 thin | 5.3041 thick 5.3141 thin | 5.3021 thick 5.3121 thin |
| 42 - 43 | 5.3002 thick 5.3102 thin | 5.3012 thick 5.3112 thin | 5.3032 thick 5.3132 thin | 5.3042 thick 5.3142 thin | 5.3022 thick 5.3122 thin |
| 44 - 45 | 5.3003 thick 5.3103 thin | 5.3013 thick 5.3113 thin | 5.3033 thick 5.3133 thin | 5.3043 thick 5.3143 thin | 5.3023 thick 5.3123 thin |
| 46 - 47 | 5.3004 thick 5.3104 thin | 5.3014 thick 5.3114 thin | 5.3034 thick 5.3134 thin | 5.3044 thick 5.3144 thin | 5.3024 thick 5.3124 thin |



blend of 80 % cotton und 20 % polyamide

hypoallergenic

washable to 60 °C

The seams are kept to a minimum and are hardly noticeable.

These socks do not have a cutting rubber edge because they are only processed with a small amount of rubber.

Measurement information for special products




Standard socks
Foot length F:
Chopart: 13 – 14 cm
Lisfranc: 19 – 20 cm



Prosthetic socks

| SIZES | 1 weiß | 2 beige | 3 marine | 4 grau | 5 schwarz |
|----------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| 35 - 38 | 5.3200 thick 5.3300 thin | 5.3210 thick 5.3310 thin | 5.3230 thick 5.3330 thin | 5.3240 thick 5.3340 thin | 5.3220 thick 5.3320 thin |
| 39 - 41 | 5.3201 thick 5.3301 thin | 5.3211 thick 5.3311 thin | 5.3231 thick 5.3331 thin | 5.3241 thick 5.3341 thin | 5.3221 thick 5.3321 thin |
| 42 - 43 | 5.3202 thick 5.3302 thin | 5.3212 thick 5.3312 thin | 5.3232 thick 5.3332 thin | 5.3242 thick 5.3342 thin | 5.3222 thick 5.3322 thin |
| 44 - 45 | 5.3203 thick 5.3303 thin | 5.3213 thick 5.3313 thin | 5.3233 thick 5.3333 thin | 5.3243 thick 5.3343 thin | 5.3223 thick 5.3323 thin |
| 46 - 47 | 5.3204 thick 5.3304 thin | 5.3214 thick 5.3314 thin | 5.3234 thick 5.3334 thin | 5.3244 thick 5.3344 thin | 5.3224 thick 5.3324 thin |



THANNER GmbH, Am Fallenweg 30, D-89420 Höchstädt/Donau, Phone +49 (0) 90 74/95 51-0, Fax +49 (0) 90 74/95 51-50
E-mail: info@thanner-gmbh.com, Internet: www.thanner-gmbh.com
Onlineshop: www.thannershop.com, [f](https://www.facebook.com/thanner.de): www.facebook.com/thanner.de